

Experienced User Protocol

Please wear gloves at all times

1. To the **PowerBead Tubes** provided, add **0.25 grams of soil sample**. (or 700ul of extracted meiofuana fraction)
2. Gently vortex to mix.
3. **Check Solution C1**. If **Solution C1** is precipitated, heat solution to 60°C until dissolved before use.
4. Add 60 µl of **Solution C1** and invert several times or vortex briefly.
5. Secure **PowerBead Tubes** horizontally using the MO BIO Vortex Adapter tube holder for the vortex (MO BIO Catalog# 13000-V1-24) or secure tubes horizontally on a flat-bed vortex pad with tape. Vortex at maximum speed for ~~10 minutes~~. **5 minutes**.
Note: If you are using the 24 place Vortex Adapter for more than 12 preps, increase the vortex time by 5-10 minutes.
6. Make sure the PowerBead Tubes rotate freely in your centrifuge without rubbing. Centrifuge tubes at 10,000 x *g* for 30 seconds at room temperature. **CAUTION:** Be sure not to exceed 10,000 x *g* or tubes may break.
7. Transfer the supernatant to a clean **2 ml Collection Tube** (provided).
Note: Expect between 400 to **500 µl of supernatant**. Supernatant may still contain some soil particles.
8. Add 250 µl of **Solution C2** and vortex for 5 seconds. Incubate at 4°C for 5 minutes.
9. Centrifuge the tubes at room temperature for 1 minute at 10,000 x *g*.
10. Avoiding the pellet, transfer up to, but no more than, 600 µl of supernatant to a clean **2 ml Collection Tube** (provided).
11. Add 200 µl of **Solution C3** and vortex briefly. Incubate at 4°C for 5 minutes.
12. Centrifuge the tubes at room temperature for 1 minute at 10,000 x *g*.
13. Avoiding the pellet, transfer up to, but no more than, 750 µl of supernatant into a clean **2 ml Collection Tube** (provided).
14. Shake to mix Solution C4 before use. Add 1200 µl of **Solution C4** to the supernatant and vortex for 5 seconds.
15. Load approximately 675 µl onto a **Spin Filter** and centrifuge at 10,000 x *g* for 1 minute at room temperature. Discard the flow through and add an additional 675 µl of supernatant to the **Spin Filter** and centrifuge at 10,000 x *g* for 1 minute at room temperature. Load the remaining supernatant onto the **Spin Filter** and centrifuge at 10,000 x *g* for 1 minute at room temperature.
Note: A total of three loads for each sample processed are required.
16. Add 500 µl of **Solution C5** and centrifuge at room temperature for 30 seconds at 10,000 x *g*.
17. Discard the flow through.
18. Centrifuge again at room temperature for 1 minute at 10,000 x *g*.
19. Carefully place spin filter in a clean **2 ml Collection Tube** (provided). Avoid splashing any **Solution C5** onto the **Spin Filter**.
20. Add 100 µl of **Solution C6** to the center of the white filter membrane. Alternatively, sterile DNA-Free PCR Grade Water may be used for elution from the silica Spin Filter membrane at this step (MO BIO Catalog# 17000-10).
21. Centrifuge at room temperature for 30 seconds at 10,000 x *g*.
22. Discard the **Spin Filter**. The DNA in the tube is now ready for any downstream application. No further steps are required.

We recommend storing DNA frozen (-20° to **-80°C**). **Solution C6** contains no EDTA. To concentrate the DNA see the Hints & Troubleshooting Guide.

Thank you for choosing the PowerSoil® DNA Isolation Kit.