

# Qubit Quantification of DNA

Eisen Lab

April 4, 2014

## 1 Introduction

The protocol uses a Qubit Fluorometer (Invitrogen) to quantify DNA, e.g. from environmental DNA extractions or purified PCR products. The Qubit kit is stored at 4°C (Quant-iT dsDNA Assay Kit, High Sensitivity, 1000 assays. Invitrogen product no. Q33120).

## 2 Protocol

### 2.1 Prepare working solution

Prepare a master mix of working solution by diluting Florescent Dye (Quant-iT dsDNA HS reagent) to Quant-iT dsDNA HS buffer at a final ratio of (1:199): add 1 µl of dye into every 198 µl of Qubit buffer. Ensure that all solutions are fully liquid before mixing; the Florescent Dye may crystallize during storage at 4°C.

### 2.2 Calibrate Qubit

Before quantifying DNA, the Qubit machine should be calibrated using standard solutions. However, if the machine has been recently calibrated, select “use last calibration” when prompted on screen.

### 2.3 Prepare DNA samples

Prepare sample solution by adding 1-10 µl of each DNA sample to working solution (190-199 µl, respectively). Each sample should be prepared in a separate 0.3 ml Qubit tube (Qubit assay tubes, Invitrogen product no. Q32856), where the final volume of the prepared sample solution is 200 µl (DNA sample mixed with working solution).

### 2.4 Measure Sample Concentration

Measure the DNA concentration of each sample, as follows:

- On the Qubit screen prompt, chose “Quant-it dsDNA, HS” (high sensitivity DNA quantification)
- Choose “use last calibration”
- Insert sample into slot, close plastic cover, and press “GO” button
- Select “calculate the sample concentration” and then choose the input sample volume used in Step 2.3
- If the Qubit cannot quantify DNA (“sample concentration too low”), prepare a new sample solution by increasing the amount of DNA diluted into the working solution (Step 2.3, for example, use 5 µl instead of 1 µl of your DNA sample). Re-measure new sample solution.